Hypro-Oss® ENGLISH Instruction for use

Hypro-Oss is a natural bovine substitution material. This patented, lyophilized, osteoinductive grafting material of animal origin type I is the result of six years of intensive research and development in cooperation between the scientific teams of the companies BIOMPLON GmbH and Hypro s.r.o.

Composition

Hypro-Oss is a natural bovine type I; it consists of approximately 70% Hydroxyapatite and 30% Atelocollagen type I.

Packaging

Hypro-Oss is packed in a sterile glass bottle and sealed in double peel-off blister packaging.

Hypro-Oss is available in the following sizes:

Grain Size 0.5 - 1 mm	Grain Size 0.5 - 2 mm	Grain Size 3 - 6 mm	Grain Size 4 - 8 mm
Cat. No. Volume	Cat. No. Volume	Cat. No. Volume	Cat. No. Volume
070 0,5 ml	074 0,5 ml	078 1 ml	084 1 ml
071 1 ml	075 1 ml	079 5 ml	085 5 ml
072 1 ml	076 3 ml	080 10 ml	086 10 ml
073 5 ml	077 5 ml	081 20 ml	087 20 ml
090 25 ml	091 25 ml	082 25 ml	088 25 ml

Hypro-Oss indication

It is used for permanent filling or reconstruction of aseptic bone defects. Indications in maxillofacial surgery and dental surgery are:

- Implantology, Periodontology, and Oral Surgery
- sinus lift
- Horizontal augmentation
- Intraosseous defects
- Post-traumatic defects
- Extraction sockets
- Vertical augmentation
- Furcation defects
- Filling of cysts
- Periodontal defects
- Filling of bone defects after excision of benign tumors and cysts.
- Indications in orthopedic surgery and traumatology as:
 - Filling of bone defects in juxta-articular changes
 - Filling of defects of the acetabulum on fracture of prosthesis
 - Filling of bone defects after excision of benign tumors
 - Filling of bone cysts
 - Filling of tissue defects in cartilage and bone transplants
 - Filling of bone defects at donor sites after harvest of autogenous bone.

Contra-indications

Acute and chronic infections at the implant site, defects in the area of open wounds, acute infectious diseases of endocrine origin, severe renal insufficiency, diseases of bone metabolism, on-going treatment with glucocorticoids and mineralocorticoids and with agents affecting calcium metabolism (e.g., calcitonin), severe or difficult to control diabetes mellitus, immunosuppressive therapy, malignancies (because of the lowered diagnostic value of the X-ray examinations in cases of tumor occurrence in the implantation site).

Administration

Hypro-Oss material is hydrophilic and it can be easily mixed with the patient's blood or sterile saline before insertion. This Hypro-Oss blood mixture has excellent coagulation properties, thus expediting extremely accurate placement. Bone regeneration can be promoted in combination with blood or sterile saline.

Side effects

After surgery intervention there may occur disturbances of wound healing and/or hematoma that may increase the risk of an infection. In cases of wound healing stability, the possibility that material may be displaced and cause local irritation cannot be excluded. There are no special instructions regarding the use of Hypro-Oss in elderly patients. It is contra-indicated to implant Hypro-Oss into the area of open epiphyseal plate in children.

Use during pregnancy and lactation

There are no studies in pregnancy or lactation available.

Precautions and warnings

Hypro-Oss is not mechanically more loaded areas, a supporting stabilization by metal osteosynthesis and plaster splint is recommended. It is recommended to verify the correct implantation with suitable imaging procedures.

Interactions

To date, there are no known interactions with drugs or other medical devices except those mentioned under Contra-indications.

Dosage and directions for use

The amount of Hypro-Oss needed depends on the size of the bone defect. Hypro-Oss is a single use device and is implanted during a surgical procedure. Hypro-Oss is intended for permanent implantation in patients.

Prerequisites for adequate implantation results

Prerequisite for an effective use is a stable Hypro-Oss implantation into a vital, well-treated bone cavity, with moderate biomechanical loading. Hypro-Oss is particularly suitable for filling cancellous bone defects and should always be grafted in a way that it will not be displaced. Hypro-Oss must be in contact with vital bone tissue. The bone defect should be completely filled with Hypro-Oss or a mixture of Hypro-Oss with autograft. The site for the implant has to be prepared by freshening the adjacent bone. The aim of a stable, intraosseous implantation is a direct osseous contact between Hypro-Oss and the vital bone without an intermediate layer of connective tissue. Hypro-Oss macro porous structure is ideal in its osteoconductive function and promotes the ingrowth of blood vessels and nerves. Native bone osteoclasts stimulates and accelerates new bone formation by triggering the thrombocytes aggregation, the degradation, and the release of growth peptides such as Beta TGF, IGF1, IGF2 and PDGF, which are responsible for the osteoconduction and consequently for the enhanced bone formation. The adhesion and spread of osteoblasts over the Hypro-Oss surface creates an open, interconnecting pore structure, thus prompting a bioactive reaction with its bone tissue formation, bone tissue strengthening and bone tissue interlinking, thus leading to the restoration of bone and its functions. Hypro-Oss is a natural bovine material, which combines the osteoinductive properties of the atelocollagen act together in harmony for enhanced, high quality new bone formation. Care shall be taken that, unless

the bone renewal process is complete, stabilization may be necessary to reduce the mechanical loading on the site of implantation (e.g., rigid dressing, metal osteosynthesis).

Hypro-Oss innovative properties

The development of the innovative Hypro-Oss product line was conceived with the idea of an ideal biomaterial in mind, a material with the highest biocompatibility and affinity to the new endogenous bone formation. Hypro-Oss is a natural bovine type I. In order to accomplish this purpose, we implemented our proprietary atelotripan and lyophilization technologies that preserve natural collagen components of bone material within bone structure after modifying the collagen to a non-immunogenic atelocollagen. This preserves the natural crystalline structure of the hydroxyapatite, as well as the collagen components by completely avoiding the heating (thermo) processing method that other manufacturers still use. Thanks to our innovative proprietary processing technology, Hypro-Oss has the following important characteristics:

- Native bovine graft components for enhanced new bone formation
- Teletropine as collagen components; non-immunogenic peptide
- Acceleration of physiological tissue healing process and bone formation
- Protects grafting site from infection (bacteriostatic effect of atelocollagen)
- Capability to apply medication to the surgical site
- Highest biocompatibility, absence of a foreign body response
- Preserves the natural structures of collagen and hydroxyapatite due to lyophilization processing.

These characteristics allow enhanced and consistent new bone formation, and persistent integration between mature newly formed bone and existing bone matrices.

BSE prevention regulations

The manufacturer issues documents at its disposal confirming that the country has established compulsory notification of BSE (Notices of the Department of Agriculture no. 287/1999 Cof. of Laws) and that all slaughtered animals older than 72 months are subject to obligatory examination based on the Notices of the Department of Agriculture no. 309/2011, Cof. of Laws. The manufacturer monitors and controls the quality of the new endogenous bone with accelerated new bone formation by the European Union, including the OIE notice Terrestrial Animal Health Code.

Deactivation or removal of infectious agents

The method of collagen treatment includes several repeated extractions that are employed to remove non-collagen globular proteins from the connective tissue, since prion (the infectious agent of BSE) belongs to the category of globular proteins. The deactivation process is an additional safeguarding procedure that results in a reduced prion content, but there is any in the tissue. Another technological process that is appropriate in view of deactivation prions is the action of a saturated calcium hydroxide solution with pH-value of 12.5. The prion protein (globular protein) contains an elongated conformation of the beta-pleated sheet type that changes to a physiologically normal conformation (non-infectious) through the action of hydroxytic substances (saturated calcium hydroxides with pH 12.5). The action of hydroxides is recommended in the WHO document CPMF and CVMF EMEA/410/01, Rev. 3 as an effective measure to reduce the risk of TSE transfer through medical devices. These procedures include structural changes such as decrystallization of prion protein (resistant to proteolysis), which results in irreversible, non-toxic conformation of the prion protein. Collagen treated this way and used in the production of Hypro-Oss consists of prion-free tissues, which are fully safe as TSE is concerned.

Sterilization method

After final packaging of Hypro-Oss's glass bottles in a double blister packaging, Hypro-Oss undergoes Gamma Irradiation as the final sterilization step. Gamma radiation sterilization is a process that effectively kills or eliminates almost all microorganisms (bacteria, viruses, and spore forms. Gamma irradiation is a physical means of decontamination because it kills bacteria by breaking down bacterial DNA, thus inhibiting bacterial division.

Shelf life

The expiry date is printed on every single packing, the primary packaging as well as the secondary packaging. Hypro-Oss must not be used after the expiration date (best before date).

Restrictions or loss of usability prior to expiration date
In case of damage to the primary packaging (bottle) or the secondary packaging (peel off packaging), the sterility of the Hypro-Oss is no longer guaranteed. Hypro-Oss must no longer be used. The content of unused but opened or damaged packages must not be re-sterilized and has to be discarded.

Storage conditions

Hypro-Oss must be stored in a dry place at room temperature. Hypro-Oss is not hygroscopic and must be stored from 0°C to +50°C. It needs to avoid sources of heat and direct sunlight.

Safe disposal

Residual Hypro-Oss is to be discarded with the hospital waste.

Classification

Hypro-Oss is a medical device class III.



Responsible worldwide distributor

Biomplon GmbH, Friedrich-List-Str. 27, 35398 Giessen, Germany, Europe
Phone +49 641-68681123, Fax +49 641-68681124, www.biomplon.de, e-mail: info@biomplon.de.

Responsible producer

Hypro s.r.o., No. 137, 763 61 Komárov, Czech Republic, E.U., phone +420 57115972, fax +420 57115974, www.hypro.cz.

Used symbols

LOT	REF	Production date	Expiry date
Batch number	Catalog number	Production date	Expiry date
			
Sterilized by gamma-radiation	For single use only	Do not sterilize	Do not use if package damaged
			
Please read the package leaflet before using	Temperature limitation	Keep away from sunlight	CE 1023 notified body, ITC Zlin, a.s., Czech Republic

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© 1.1.2015 Hypro s.r.o.
Biomplon GmbH, Friedrich-List-Str. 27, 35398 Giessen, Bundesrepublik Deutsch