kasios®

kasios® TCP

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Synthetic phosphocalcic ceramic made from β tri-calcium phosphate (β TCP) which closely resembles the mineral phase of natural bone.

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Kasios®TCP is entirely synthetic, so that risks as regards immune response or infection are eliminated.

Cavity filling and graft integration radiological monitoring is not impeded since the TCP granules are radio-opaque.

The shape of Kasios®TCP granules solves the problem of filling irregularly shaped cavities. Because of its composition, Kasios®TCP is biocompatible. Osteo-integration and resorption proceed at a rapid pace.

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Réf. / Ref.	Description / Description	Conditionnement / Packaging
K4070100	Granules 150-500µm	Boîte de 2 doses stériles de 1g / Box of 2 sterile 1g dose
K4070150 ✓	Granules 150-500µm	Boîte de 5 doses stériles de 1g / Box of 5 sterile 1g dose
K4070500	Granules 500-1000µm	Boîte de 5 doses stériles de 1g / Box of 5 sterile 1g dose
K4071000	Granules 1000-2000µm	Boîte de 5 doses stériles de 2g / Box of 5 sterile 2g dose

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Hold the blister firmly in one hand and slowly peel off the seal, being careful to hold at a slant to avoid spilling granules.

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
Mix the Kasios[®] TCP granules with either the patient's blood, saline solution, or distilled water. Do not humidify the granules excessively. Any surplus liquid should be wiped away with the corner of a sterile sponge.

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The flat side of a sterile spatula may be used to mix Kasios[®] TCP granules and the liquid. Humidified granules can be transported to the surgical area using the curved part of the spatula.

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Gently fill each defect progressively up to the highest level of the bone cavity. Osteo-integration will be improved if blood and vessels can penetrate the graft.

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Once filling with the Kasios[®] TCP granules is completed, the flap should be carefully adjusted and completely closed. A membrane can also be used to ensure full closure.

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